

## Blood and Body Fluid DNA Extraction for Streptococci

(Carvalho Mda G et al. Evaluation and improvement of real-time PCR assays targeting *lytA*, *ply*, and *psaA* genes for detection of pneumococcal DNA. J Clin Microbiol. 2007 Aug;45(8):2460-6. Epub 2007 May 30.)

- *When using a new QIAampDNA Mini kit or DNeasyBlood & Tissue Kit, add alcohol to buffers AW1 and AW2 as indicated on the bottles.*
- *Equilibrate samples and AE Buffer to room temperature (15-25 °C).*
- *Heat the water bath to 56 °C.*
- *If a precipitate has formed in Buffer AL, dissolve by incubating at 56 °C.*
- *All centrifugations steps should be carried out at room temperature.*
  
- Pipet **50 uL** of TE buffer containing **0.08g/mL** of lysozyme + **150 U/mL** of mutanolysin into the bottom of a 1.5 mL microcentrifuge tube.
  - **TE buffer (10mM Tris-HCl, 1 mM EDTA, pH 8.0)**
  - **The lysozyme (Sigma-L-6876), and mutanolysin (Sigma-M9901) should be added to the TE buffer right before the use (freshly prepared).**
- Add 200 uL of the sample. Vortex and incubate at 37 °C in water bath for 1 hour.
- Add 20 uL of proteinase K. Vortex tube briefly.
- Incubate at 56 °C in water bath for 30 min.
- Centrifuge briefly to remove drops from inside of the lid.
- Add 200 uL of Buffer AL. Vortex tube briefly.
- Incubate at room temperature for 10 min.
- **Add 250 uL ethanol (96-100%) to the sample and mix gently pipetting up and down avoiding splash.**
- Pipet the mixture into the QIAamp Spin column sitting in the 2 mL collection tube without wetting the rim, close the cap.
- Centrifuge at 6000 x g or 8000 rpm for 1 min. Discard eluate and collection tube.

- Place mini column into fresh collection tube. Add 500  $\mu$ L of Buffer AW1 without wetting the rim, close the cap and centrifuge 1 min. at 6000 x  $g$  or 8000 rpm. Discard eluate and collection tube.
- Place mini column into fresh collection tube. Add 500  $\mu$ L Buffer AW2 and centrifuge for 3 minutes at FULL SPEED (20,000 x  $g$  or 14,000 rpm) to dry column. Discard eluate and collection tube.
- Place mini column in 1.5 ml sterile eppendorf tube. Add 100  $\mu$ L Buffer AE onto membrane. Incubate room temperature 5 minutes. Centrifuge 1 min. at 6000 x  $g$  or 8000 rpm to eluate DNA. (Keep at -20  $^{\circ}$ C until use).